

Biotoools

Mechanical Pipettes

improve accuracy and reproducibility of your experiments

Biotoools offers a portfolio of **Mechanical Pipettes** of Adjustable Volumes.

Biotoools single-channel pipettes are available in eight volume ranges for reliable dispensing of volumes from 0.1 μ l to 1 ml.

- n mechanical pipette with adjustable volume
- n high precision light
- n weight ergonomic design
- n excellent quality-price ratio
- n easy to calibrate and maintain
- n supplied with a safety filter

Ordering information

Cat. No.Format

50.201	Pipette B 10: 0.5-10 μ l, increment 0.1 μ l
50.202	Pipette B 20: 2-20 μ l, increment 0.5 μ l
50.203	Pipette B 50: 5-50 μ l, increment 0.5 μ l
50.204	Pipette B 100: 10-100 μ l, increment 1 μ l
50.205	Pipette B 200: 20-200 μ l, increment 1 μ l
50.206	Pipette B 200: 50-200 μ l, increment 1 μ l
50.207	Pipette B 1000: 100-1000 μ l, increment 5 μ l
50.208	Pipette B 1000: 200-1000 μ l, increment 5 μ l



Additional information at info@biotoools.eu

Biotoools B&M Labs, S.A.
Valle de Tobalina, 52 Nave 39
28021 Madrid, Spain
Tel.: +34 917 100 074
Fax: +34 915 053 118
info@biotoools.eu



BIOTOOLS
BIOTOOLS B&M LABS.S.A.