

## LionXT

	LionXT 75	LionXT 100	LionXT 150
<b>Lenght</b>	76 cm	106 cm	160 cm
<b>Depth</b>	52 cm	70 cm	70 cm
<b>Height</b>	66 cm	72 cm	72 cm
<b>Capacity</b>	12 microplates	20 microplates	36 microplates
<b>Arms</b>	1	1 or 2	1 or 2
<b>Channels</b>	1 or 2	2, 4, 8 or 96	2, 4, 8 or 96
<b>Pump</b>	2-250 $\mu$ 10-1000 $\mu$	2-250 $\mu$ 10-1000 $\mu$	2-250 $\mu$ 10-1000 $\mu$
<b>Dispense</b>	Independent each channel	Independent each channel	Independent each channel
<b>Liquid detection</b>	Pressure	Pressure	Pressure
<b>Data</b>	Import and export data from excel		
<b>Main characteristics</b>	Independent movement in X and Y axis		
<b>Extraction method</b>	Magnetic Bead Station -	Magnetic Bead Station SPE Vacuum Station	Magnetic Bead Station SPE Vacuum Station
<b>Type</b>	DNA RNA Protein	DNA RNA Protein	DNA RNA Protein
<b>Optional Devices</b>	HEPA filter Plate Handler "Plate-GripX" Barcode Scanner Plate Shaker "Plate-ShakeX" Heating/Cooling Devices	HEPA filter Plate Handler "Plate-GripX" Barcode Scanner Plate Shaker "Plate-ShakeX" Heating/Cooling Devices	HEPA filter Plate Handler "Plate-GripX" Barcode Scanner Plate Shaker "Plate-ShakeX" Heating/Cooling Devices