



Oligo dT(15)

Product Description:

A primer designed to initiate synthesis of a cDNA from total RNA in a reverse transcription reaction. For use in generating labeled cDNA for screening microarrays.

Form: Supplied in lyophilized form

Storage: -20°C

T_m: 26°C

Purity:

Purified by reverse phase chromatography, >95%

Pack size:

30 µg equals ~ 1 OD unit (A260)

Quality Control:

Purity is verified by analytical PAGE

Notes for Use

For maximum recovery prior removing the cap of the tube spin it briefly to bring down the content. Reconstitute with DNase free water and vortex briefly. Avoid repeat freeze-thaw cycles.

Order Information:

Cat. No.	Product	
20.401	Random Primers	30 µg
20.402	Oligo dT(15)	30 µg

Notice to buyers/users:

Some of the applications which may be performed with this product are covered by applicable patents in certain countries. The purchase of this product does not include or provide a license to perform patented applications. Users may be required to obtain a license depending on the country and/or application.



Oligo dT(15)

Descripción del Producto:

Primer diseñado para iniciar la síntesis de ADNc a partir de ARN total en reacciones de transcripción reversa. Se puede emplear en marcaje general de ADNc en screening de microarrays.

Formato: liofilizado

Almacenamiento: -20°C

MW: 4501 mol/l

Tm: 26°C

Pureza:

Purificado mediante cromatografía en fase reversa, >95%

Tamaño del Producto:

30 µg ~ 1 OD unit (A260)

Control de Calidad:

La pureza del producto se ha verificado mediante analitical PAGE

Notas de Uso

Spin el vial brevemente antes de abrirlo para así recuperar todo su contenido. Reconstituya el vial con agua libre de DNAsas y mezcle con el vortex brevemente. Evite ciclos de congelación y descongelación.

Información de Pedidos:

Cat. No.	Producto	
20.401	Random Primers	30 µg
20.402	Oligo dT(15)	30 µg

Nota a compradores/usuarios:

Algunas de las aplicaciones que pueden ser realizadas con este producto están protegidas por patentes aplicables en ciertos países. La compra de este producto no incluye una licencia para llevar a cabo aplicaciones patentadas. Los usuarios podrían tener que adquirir una licencia dependiendo del país y/o aplicación.